



National Pharmaceutical Control Bureau  
MINISTRY OF HEALTH MALAYSIA



WHO Collaborating Centre  
for Regulatory Control of  
Pharmaceuticals



Pharmaceutical Inspection  
Convention and Pharmaceutical  
Inspection Co-operation  
Scheme



SIRIM  
Certified to ISO 9001:2000  
Cert. No: AR 2293



MS ISO/IEC 17025:2005  
NO. SAKM 450

# MICROBIAL CONTAMINATION TEST (MCT)

---

Centre for Quality Control

National Pharmaceutical Control Bureau

Lot 36, Jalan Universiti, 46200 Petaling Jaya, Selangor

DL: +6.03.78835400 (EXT5442) | F: +6.03.79567075 |

WS : [www.bpfk.gov.my](http://www.bpfk.gov.my) |

---



# OUTLINE

- **Introduction**
- **Certificate of Analysis**
- **Media Validation**
- **Test Method**
  - Total Viable Aerobic Count
  - Test for Specified Microorganisms
- **Method Validation**



# Introduction - Microbial Contamination Test (MCT)

- **Microbial Contamination Test is conducted on non-sterile products to check:**
  - The level of microbial (bacterial and fungal) contamination
  - Presence/ absence of certain pathogenic microorganism **in order to assure product safety.**
- **Types of samples include:**



Capsule



Tablet



Aqueous preparation



Transdermal Patch



Cream



Pessary



Inhaler



Suppository



# Certificate of Analysis

## ❑ **Specification and results**

- refer **British Pharmacopoeia 2012**,  
Table 5.1.4-1 Acceptance criteria for microbiological  
quality of non-sterile dosage forms




Table 5.1.4-1. – Acceptance criteria for microbiological quality of non-sterile dosage forms

Route of administration	TAMC (CFU/g or CFU/mL)	TYMC (CFU/g or CFU/mL)	Specified micro-organisms
Non-aqueous preparations for oral use	10 <sup>3</sup>	10 <sup>2</sup>	Absence of <i>Escherichia coli</i> (1 g or 1 mL)
Aqueous preparations for oral use	10 <sup>2</sup>	10 <sup>1</sup>	Absence of <i>Escherichia coli</i> (1 g or 1 mL)
Rectal use	10 <sup>3</sup>	10 <sup>2</sup>	-
Oromucosal use Gingival use Cutaneous use Nasal use Auricular use	10 <sup>2</sup>	10 <sup>1</sup>	Absence of <i>Staphylococcus aureus</i> (1 g or 1 mL) Absence of <i>Pseudomonas aeruginosa</i> (1 g or 1 mL)
Vaginal use	10 <sup>2</sup>	10 <sup>1</sup>	Absence of <i>Pseudomonas aeruginosa</i> (1 g or 1 mL) Absence of <i>Staphylococcus aureus</i> (1 g or 1 mL) Absence of <i>Candida albicans</i> (1 g or 1 mL)
Transdermal patches (limits for one patch including adhesive layer and backing)	10 <sup>2</sup>	10 <sup>1</sup>	Absence of <i>Staphylococcus aureus</i> (1 patch) Absence of <i>Pseudomonas aeruginosa</i> (1 patch)
Inhalation use (special requirements apply to liquid preparations for nebulisation)	10 <sup>2</sup>	10 <sup>1</sup>	Absence of <i>Staphylococcus aureus</i> (1 g or 1 mL) Absence of <i>Pseudomonas aeruginosa</i> (1 g or 1 mL) Absence of bile-tolerant gram-negative bacteria (1 g or 1 mL)
♦Special Ph. Eur. provision for oral dosage forms containing raw materials of natural (animal, vegetal or mineral) origin for which antimicrobial pretreatment is not feasible and for which the competent authority accepts TAMC of the raw material exceeding 10 <sup>3</sup> CFU/g or CFU/mL.	10 <sup>4</sup>	10 <sup>2</sup>	Not more than 10 <sup>2</sup> CFU of bile-tolerant gram-negative bacteria (1 g or 1 mL) Absence of <i>Salmonella</i> (10 g or 10 mL) Absence of <i>Escherichia coli</i> (1 g or 1 mL) Absence of <i>Staphylococcus aureus</i> (1 g or 1 mL)♦



# Media Validation



## Prior to test, make sure that:

- ✓ Media is sterile
- ✓ Media supports growth of microorganisms
- ✓ Selective media is selective  
(promote certain organisms but inhibit non-target organisms)


## In order to so,

- Test for Media Sterility
- Test for Growth Promotion & Inhibitory Properties



# Media Validation- Test for Media Sterility

- To prevent False Positive result**
  - maybe due to contaminated media
- To ensure the media is sterile**
- Negative Control**
  - Use the chosen sterile diluents in place of the sample under test
  - Alternatively, incubate portions of the media for a few days at the specified temperature.
- Acceptance criteria:** No growth observed





# Media Validation- Test for Growth Promotion and Inhibitory Properties

□ **There are 2 categories of media used in MCT:**

## **1. General nutritive media**

- used in Total Viable Aerobic Count
- suitable for cultivation of a wide variety of microorganisms
- e.g. Tryptone Soya Agar
- Test for Growth Promotion Properties

## **2. Selective media**

- used in Test for Specified Microorganisms
- contains ingredients which promotes growth of certain organisms but inhibit other non-target microorganisms
- e.g. Mannitol Salt Agar, Cefrimide Agar, MacConkey Broth
- Test for Growth Promotion, Indicative and Inhibitory Properties






# Media Validation- General Nutritive Media

## □ Test for Growth Promotion Properties

- ✓ To verify that media used are able to support growth of a wide variety of microorganisms

### Test Method

- Inoculate portions/ plates of media with a small number ( $< 100$  cfu) of microorganisms\* indicated in Table 1.
- Use a separate plate of medium for each microorganism.
- Incubate at the specified temperature.



\*Note: Microorganisms used should not be more than 5 passages removed from the original seed-lot.



# Media Validation- General Nutritive Media

**Table 1- Media, Microorganisms and Test Condition for Growth Promotion Test**

Test	Media Used	Microorganisms	Test Condition
<b>Total Aerobic Microbial Count (TAMC)</b>	Tryptone Soya Agar (TSA)	<ul style="list-style-type: none"><li>• Staphylococcus aureus</li><li>• Pseudomonas aeruginosa</li><li>• Bacillus subtilis</li><li>• Candida albicans</li><li>• Aspergillus brasiliensis</li></ul>	$\leq 100$ cfu 30 - 35°C, $\leq 3$ days for bacteria and $\leq 5$ days for fungi
	Tryptone Soya Broth (TSB)	<ul style="list-style-type: none"><li>• Staphylococcus aureus</li><li>• Pseudomonas aeruginosa</li><li>• Bacillus subtilis</li></ul>	$\leq 100$ cfu 30 - 35°C, $\leq 3$ days
<b>Total Yeasts and Moulds Count (TYMC)</b>	Sabouraud Dextrose Agar (SDA)	<ul style="list-style-type: none"><li>• Candida albicans</li><li>• Aspergillus brasiliensis</li></ul>	$\leq 100$ cfu 20 - 25°C, $\leq 5$ days



# Media Validation- General Nutritive Media


## Acceptance Criteria

### **Solid media:**

- Growth obtained must not differ by a factor of 2 (50-200%) from the calculated value for a standardized inoculum. (Quantitative)
- Growth of the microorganisms comparable to that previously obtained with a previously tested and approved batch of medium occurs.

### **Liquid media:**

- Clearly visible growth of microorganisms comparable to that previously obtained with a previously tested and approved batch of medium occurs



# Media Validation - Selective Media


- For media used in Test for Specified Microorganisms
- Tests for Growth Promotion, Indicative and Inhibitory Properties need to be conducted

## 1. Test for Growth Promoting Properties

### Liquid & Solid Media

1. Inoculate a portion of the medium with a small number ( $\leq 100$  cfu) of the appropriate microorganism (Table 2). For Solid media, use surface spread method.
2. Incubate at the specified temperature for not more than the shortest time specified in the test.

Acceptance criteria: Clearly visible growth



# Media Validation - Selective Media

## 2. Test for Inhibitory Properties

1. Inoculate the medium with at least 100 cfu of the appropriate microorganism (Table 2).
2. Incubate at the specified temperature for not less than the longest time specified in the test.

Acceptance criteria: No Growth of the test microorganisms occurs

## 3. Test for Indicative Properties

1. Inoculate each plate of medium using surface spread method with a small number ( $\leq 100$  cfu) of the appropriate microorganism (Table 2).
2. Incubate at the specified temperature for a period of time within the range specified in the test.

Acceptance criteria: Colonies are comparable in appearance and indicative reactions to those previously obtained with a previously tested and approved batch of medium.



# Media Validation - Selective Media

**Table 2- Growth Promoting, Inhibitory and Indicative Properties of Media**

Test for	Media	Property	Test Strain
Bile-Tolerant Gram Negative Bacteria	Enterobacteria Enrichment Broth (EEB)	Growth Promoting	E. coli P. aeruginosa
		Inhibitory	S. aureus
	Violet Red Bile Glucose Agar (VRBGA)	Growth Promoting & Indicative	E. coli P. aeruginosa
Escherichia coli	MacConkey Broth (MCB)	Growth Promoting	E. coli
		Inhibitory	S. aureus
	MacConkey Agar (MCA)	Growth Promoting & Indicative	E. coli
Salmonella	Rappaport Vassiliadis Salmonella Enrichment Broth (RVS)	Growth Promoting	Salmonella typhimurium or Salmonella abony
		Inhibitory	S. aureus
	Xylose, Lysine Deoxycholate Agar (XLD)	Growth Promoting & Indicative	Salmonella typhimurium or Salmonella abony
Pseudomonas aeruginosa	Cetrimide Agar (CETA)	Growth Promoting	P. aeruginosa
		Inhibitory	E. coli
Staphylococcus aureus	Mannitol Salt Agar (MSA)	Growth Promoting & Indicative	S. aureus
		Inhibitory	E. coli
Candida albicans	Sabouraud Dextrose Broth (SDB)	Growth Promoting	C. albicans
	Sabouraud Dextrose Agar (SDA)	Growth Promoting & Indicative	C. albicans



# Test Method

## ☐ **MCT consists of 2 tests:**

### **1. Total Viable Aerobic Count (TVAC)**

- Enumeration of bacteria and fungi present in the product
- Total Aerobic Microbial Count (TAMC)
- Total Yeast and Mould Count (TYMC)

### **2. Test for Specified Microorganism**

- Qualitative: Presence or absence of specified microorganisms
- Semi Quantitative: Test for Bile-Tolerant Gram Negative Bacteria

**\*The type of specified microorganisms tested depends of the route of administration and the type of preparation**



# Test Method - Total Viable Aerobic Count (TVAC)

- **The choice of method is based on factors such as the nature of product and the required limit of microorganisms.**

Membrane Filtration	Plate Count	Most Probable Number (MPN)
Suitable for soluble and filterable samples	Surface Spread & Pour Plate	Low precision and accuracy
Filter pore size $\leq 0.45 \mu\text{m}$	Perform test at least in duplicate for each medium	Only for Total Aerobic Microbial Count (TAMC)
Bacteria retaining efficiency of filter not affected by sample	Take arithmetic mean count for each medium	May be suitable for samples with very low bioburden





# Test Method- TVAC Membrane Filtration



- Use sterilized filtration apparatus.
- Membrane pore size  $\leq 0.45\mu\text{m}$ .




- Filter sample preparation containing 1g of product.
- Rinse the filter with an appropriate volume of diluent.



- Transfer the membrane filter to the surface of TSA and SDA for enumeration of TAMC and TYMC respectively.
- Incubate TSA at 30 - 35°C for  $\geq 3$  days and SDA at 20 - 25°C for  $\geq 5$  days.



# Test Method - TVAC Plate Count Method






# Test Method - Test for Specified Microorganisms

**Specified microorganisms tested for in MCT are...**

**Escherichia coli**



**Pseudomonas  
aeruginosa**



**Staphylococcus  
aureus**



**Salmonella**



**Candida albicans**



# Test Method - Test for *Staphylococcus aureus* & *Pseudomonas aeruginosa*

Preparations tested include those of orocumosal, gingival, cutaneous, nasal, auricular and vaginal use and transdermal patches.

10 g sample



90ml of Buffered NaCl Peptone Solution

10 ml



90ml of Tryptone Soya Broth (TSB),  
30 -35°C, 18 – 24hrs

Subculture on



Mannitol Salt Agar (MSA)

30 -35°C, 18 – 72hrs



Cetrimide Agar (CETA)



# Test Method - Test for Escherichia coli

Preparations tested include aqueous and non- aqueous preparations, oral dosage forms containing natural origin and solely herbal medicinal products.

10g/10 ml sample



90ml of Buffered NaCl Peptone Solution

10 ml



90ml of Tryptone Soya Broth (TSB),  
30 -35°C, 18 – 24hrs

1 ml



100ml of MacConkey Broth (MCB),  
42 - 44°C, 18 – 72hrs



Subculture on




30 -35°C, 18 – 72hrs

MacConkey Agar (MCA)



# Test Method - Test for *Candida albicans*

Preparation tested  
is those of vaginal  
use.



**10 g sample**



**90ml of Buffered NaCl  
Peptone Solution**

**10 ml**



**100ml of Sabouraud Dextrose Broth (SDB),  
30 -35°C, 3 – 5 days**

**Subculture on**



**30 -35°C, 24- 48hrs**

**Sabouraud Dextrose Agar (SDA)**



NPCB  
MOH



# Test Method - Test for Bile Tolerant Gram-Negative Bacteria

Preparations tested include products for inhalation, oral dosage forms containing natural origin and solely herbal medicinal products.

10 g/ 10ml sample






90ml of Tryptone Soya Broth (TSB),  
20 -25°C, 2 - 5hrs



**Enterobacteria Enrichment  
Broth- Mossel (EEB),  
30 -35°C, 24 - 48hrs**

Subculture on

30 -35°C, 18 - 24hrs




**Violet Red Bile Glucose Agar (VRBGA)**



# Test Method - Test for Specified Microorganism

- If growth is observed on selective agar, gram stain and identify the bacteria







# Method Validation

- Also known as 'Suitability of the Counting Method/ Test Method'
- To establish the ability of the chosen test method to detect microorganisms in the presence of product
- Product specific, i.e. need to conduct MCT validation on every product
- If the product contains antimicrobial ingredient/activity (e.g. antibiotic, preservative), this should be insofar possible removed or neutralised
- If surface active substance are used for sample preparation, their absence of toxicity for organisms and their compatibility with inactivators must be demonstrated
- Suitability must be confirmed if any changes which may affect the test outcome is introduced (e.g. change in formulation, change in API or preservative content)

## How?

Spike a small number of microorganisms into the product, run the test as per the chosen method, and check if the method is able to recover the microorganisms



# Validation of Total Viable Aerobic Count by Plate Count Method

**Objective:** To demonstrate the ability of the test method to detect microorganisms present in the product

- Conducted in the presence & absence of product
- Spiked known number of microorganisms  
(to obtain an inoculum of **NMT 100 CFU**. The volume of the suspension of the inoculum **should not exceed 1%** of the volume of diluted product.)

$$\text{Recovery} = \frac{\text{mean no. of colonies in presence of product}}{\text{mean no. of colonies in absence of product}} \times 100\%$$

**Acceptance Criteria:** Mean count of any test organisms not differing by a factor greater than 2 (50% – 200%)



# Validation of Total Viable Aerobic Count by Plate Count Method

In presence  
of product



10 g/ 10ml sample



90ml buffered NaCl- peptone (1:  
10)  
**+ 1 ml of microorganisms  
suspension**



TSA / SDA  
(duplicate)

Incubate

Count the mean no. of  
colonies on plates

10ml diluents

In absence of  
product



90ml buffered NaCl- peptone (1:  
10)  
**+ 1 ml of microorganisms  
suspension**



TSA / SDA  
(duplicate)

Incubate

Count the mean no. of  
colonies on plates



# Validation of Total Viable Aerobic Count by Plate Count Method

In presence  
of product



10 g/ 10ml sample



90ml buffered NaCl-peptone (1: 10)

1 ml

1 ml



~~TSA / SDA (duplicate)  
+ 10 µl microorganisms  
suspension~~

Incubate

Count the mean no. of  
colonies on plates




# Validation of Total Viable Aerobic Count by Plate Count Method

In presence  
of product



10 g/ 10ml sample




90ml buffered NaCl-peptone (1:  
10)



10 ml



+ 100  $\mu$ L microorganisms  
suspension



TSA / SDA (duplicate)

↓ Incubate

Count the mean no. of  
colonies on plates





# Validation of Total Viable Aerobic Count by Plate Count Method

## British Pharmacopoeia 2012:

### *4-5 Suitability of the counting method in the presence of product*

**4-5-1 Preparation of the sample** The method for sample preparation depends upon the physical characteristics of the product to be tested. If none of the procedures described below can be demonstrated to be satisfactory, an alternative procedure must be developed.



**4-5-2 Inoculation and dilution** Add to the sample prepared as described above (4-5-1) and to a control (with no test material included) a sufficient volume of the microbial suspension to obtain an inoculum of not more than 100 CFU. The volume of the suspension of the inoculum should not exceed 1 per cent of the volume of diluted product.

To demonstrate acceptable microbial recovery from the product, the lowest possible dilution factor of the prepared sample must be used for the test. Where this is not possible due to antimicrobial activity or poor solubility, further appropriate protocols must be developed. If inhibition of growth by the sample cannot otherwise be avoided, the aliquot of the microbial suspension may be added after neutralisation, dilution or filtration.



**4-5-4-2 Plate-count methods** Perform plate-count methods at least in duplicate for each medium and use the mean count of the result.

#### **4-5-4-2-1 Pour-plate method**

For Petri dishes 9 cm in diameter, add to the dish 1 mL of the sample prepared as described under 4-5-1 to 4-5-3 and 15-20 mL of casein soya bean digest agar or Sabouraud-dextrose agar, both media being at not more than 45 °C. If larger Petri dishes are used, the amount of agar medium is increased accordingly. For each of the micro-organisms listed in Table 2.6.12.-1, at least 2 Petri dishes are used. Incubate the plates as indicated in Table 2.6.12.-1. Take the arithmetic mean of the counts per medium and calculate the number of CFU in the original inoculum.



# Validation of Total Viable Aerobic Count by Plate Count Method

## Media, Microorganisms and Test Condition for Validation of Total Viable Aerobic Count

Test	Media Used	Microorganisms	Test Condition
<b>Total Aerobic Microbial Count (TAMC)</b>	Tryptone Soya Agar (TSA)	<ul style="list-style-type: none"><li>• Staphylococcus aureus</li><li>• Pseudomonas aeruginosa</li><li>• Bacillus subtilis</li><li>• Candida albicans</li><li>• Aspergillus niger</li></ul>	$\leq 100$ cfu 30 - 35°C, $\leq 3$ days for bacteria and $\leq 5$ days for fungi
<b>Total Yeasts and Moulds Count (TYMC)</b>	Sabouraud Dextrose Agar (SDA)	<ul style="list-style-type: none"><li>• Candida albicans</li><li>• Aspergillus niger</li></ul>	$\leq 100$ cfu 20 - 25°C, $\leq 5$ days



# Validation for Test for Specified Microorganisms

## Example: Test for Specified Microorganisms in Topical Preparation

The specified microorganisms tested are:

- i) *Staphylococcus aureus*
- ii) *Pseudomonas aeruginosa*

### Acceptance criteria:

*S. aureus* & *Ps. Aeruginosa* must be detected



10 g/ 10ml sample



90ml of Buffered NaCl Peptone Solution

+  $\leq 100$  cfu *Ps. aeruginosa*

10 ml



**Enrichment:**

90ml of Tryptone Soya Broth (TSB), 30 - 35°C, 18 - 24hrs

Subculture on



30 -35°C, 18 - 72hrs

**Selective Agar: Cetrinide Agar (CETA)**





# Validation for Test for Specified Microorganisms


## Example: Test for Specified Microorganisms in Pharmaceutical Oral Dosage Forms

The specified microorganisms tested are:

- i) Escherichia coli

### Acceptance Criteria:

E. coli must be detected





# Checklist for MCT

Test	Document required	Method	Results (Raw data)
CoA	1. Specification and Results	-	-
Routine Test	1. Total Viable Aerobic Count (TAMC and TYMC)	√	√
	2. Test for Specified Microorganism	√	√
	3. Test for Growth Promoting, Indicative and Inhibitory Properties of Media	√	√
	4. Test for Media Sterility	√	√
Validation Test	1. Total Viable Aerobic Count (TAMC and TYMC)	√	√
	2. Test for Specified Microorganism	√	√



# Comments for MCT (BM)

## Ujian Kontaminasi Mikrobial (MCT):

1. Sila kemukakan tatacara pengujian (SOP) dan keputusan ujian (raw data) untuk yang berikut:
  - Test for Growth Promoting and Inhibitory Properties dan Media Sterility Test bagi semua media yang digunakan.
  - Total Viable Aerobic Count (TAMC & TYMC)
  - Test for Specified Microorganisms

Tatacara hendaklah spesifik kepada produk. Salinan terus dari farmakopoeia tidak diterima.

2. Sila kemukakan tatacara validasi untuk ujian Total Viable Aerobic Count & Test for Specified Microorganisms, berserta acceptance criteria dan keputusan dalam bentuk raw data yang menunjukkan bahawa kandungan produk ini tidak merencatkan pertumbuhan mikroorganisma semasa MCT dijalankan.

(Sila rujuk British Pharmacopoeia – Suitability of the Counting Method in the Presence of Product & Suitability of the Test Method)

Kesemua raw data yang dikemukakan perlu mengandungi nama dan nombor kelompok bagi Finished Product, tarikh mula dan selesai pengujian, keputusan pemerhatian setiap hari & tandatangan/ nama penganalisis.

3. Sila kemukakan terjemahan bahasa Inggeris sekiranya data adalah dalam bahasa negara asing.



# Comments for MCT (English)

## Microbial Contamination Test (MCT):

1. Please provide method (SOP) and result in raw data for below:
  - Test for Growth Promoting and Inhibitory Properties dan Media Sterility Test for all the media used.
  - Total Viable Aerobic Count (TAMC & TYMC)
  - Test for Specified Microorganisms

Method must be specific to the product and photocopy from pharmacopoeia is not acceptable.

2. Please provide the validation method for Total Viable Aerobic Count & Test for Specified Microorganisms, together with acceptance criteria and the result in raw data.

(Please refer to British Pharmacopoeia – Suitability of the Counting Method in the Presence of Product & Suitability of the Test Method)

All the raw data provided must include product's name, batch number of finished product, starting date and finishing date, observation result in interval period, analyst's name and signature.

3. Please translate into English or BM if the raw data provided are in others language.



NPCB  
MOH

**THANK YOU!**